

# STROBOSCOPE

Product Codes: 08235

*Compact, battery-operated stroboscope enables operator to view devices while they are cycling.*

The rechargeable handheld device provides mobility throughout the lab. It delivers fast and accurate cycle per minute measurements to verify cycle rates.

- User adjustable flash rates - 40.0 to 12,500 fpm (flashes per minute)
- Advance or retard flash timing by degrees of angle, without changing actual fpm
- Includes AC adapter and extra flash bulb



No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form or by any means, electronic, mechanical, photocopying, recording, scanning or otherwise without the prior written permission of the Publisher.

**Limit of Liability/Disclaimer of Warranty:** While the publisher and author have used their best efforts in preparing this document, they make no representation or warranties with respect to the accuracy or completeness of the contents of this document and specifically disclaim any implied warranties of merchantability or fitness for a particular purpose. No warranty may be created or extended by sales representatives or written sales material. The advice and strategies contained herein may not be suitable for your situation. Neither the Publisher nor author shall be liable for any loss of profit or any other commercial damages, including but not limited to special, incidental, consequential or other damages.

455 Boleskine Road, Victoria B.C. Canada, V8Z 1E7  
info@ViVitroLabs.com | (250) 388-3531 | [www.ViVitroLabs.com](http://www.ViVitroLabs.com)

Copyright © 2016 by ViVitro Labs Inc. All rights reserved.

**ViVitro Labs Inc.**  
First in Cardiovascular Device Testing