

HITEST SOFTWARE


Product Code: 05831

ViVitro offers a family of quality data acquisition software products to aid in collection of data to meet ISO5840 requirements including ViVitest, Hitest and QCtest. The software is pre-installed on a laptop and tested with all hardware and A-D converters prior to shipment. Each complete system is supplied with the following components to make up a functional system:

- I/O Module – inputs/outputs signals from the software
- Connecting cables
- Software (ViVitest, Hitest or QCtest)
- AmPack Pressure Measuring System including transducers and amplifier

When coupled with the HiCycle or the Real-Time Wear Tester, Hitest provides the most powerful, user-friendly tool for durability testing. Operators are able to monitor pressure and collect data through an intuitive user interface. Hitest provides accurate pressure data and statistical analysis to ensure the tested devices are meeting ISO5840 standards.

The results highlight the percent of cycles and percent time per cycle the device is meeting or exceeding target pressure. Hitest has been designed and validated to meet the stringent requirements of ISO 5840.



CONTACT US

HEAD OFFICE

ViVidro Labs Inc.
455 Boleskine Road
Victoria, BC, V8Z 1E7, Canada
1 250.388.3531
info@vivitrolabs.com
skype: Vivitro.labs.inc.
www.vivitrolabs.com

NORTH AMERICA

Douglas Bigrigg, B. Eng.
Major Accounts Manager

ViVidro Labs Inc.
455 Boleskine Road
Victoria, BC, V8Z 1E7, Canada
1 250.388.3531 x217
DBigrigg@vivitrolabs.com
www.vivitrolabs.com

LATIN AMERICA (EXCLUDING MEXICO)

Andreluiz de Sousa Cunha

Crions

R. Viamão 1267,
LJ8 Belo Horizonte, Brazil
30431-253 Brazil
+55(31) 2515-2505 / 2515-7017
andre@crions.com
www.crions.com

EUROPE, ISRAEL AND AFRICA

Karim Mouneimne
President & CEO

ProtomedLabs SASU

24-26 BD Gay Lussac
Immeuble Energie
F-13014 Marseille, France
+33 4 86 68 68 10
k.mouneimne@protomedlabs.com
www.protomedlabs.com

CHINA

Shining Y. W. Sun
沈阳福帝斯科技有限公司

Shenyang Fortis Tech. Ltd.

024-23286966
yuwei.sun@fortis-hitech.com
www.fortis-hitech.com

JAPAN

Toshiharu Miyake
CEO

Comfort Medical, Inc

Omagi 1751-1-101
Midori-ku, Saitama-shi
Saitama, Japan
+33 4 86 68 68 10
miyake@comfort-md.com
www.comfort-md.com

MIDDLE EAST

Mohan Das

Middle East Medical Supplies

PO Box 26143 Tahir Musbah Al Kendi Bldg.
Abu Dhabi, United Arab Emirates
+971 506221598
mohanda@memsuae.com
www.memsuae.com

INDIA

Harish Rao

Emergent Medi-Tech (I) Pvt. Ltd.

116-119 Damji Shamji Industrial Complex
9, L.B.S. Marg, Kurla (West)
Mumbai, 400070, India
91-98-20210150
91-22-67825555
info@emergent.in
www.emergent.in

No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form or by any means, electronic, mechanical, photocopying, recording, scanning or otherwise without the prior written permission of the Publisher.

Limit of Liability/Disclaimer of Warranty: While the publisher and author have used their best efforts in preparing this document, they make no representation or warranties with respect to the accuracy or completeness of the contents of this document and specifically disclaim any implied warranties of merchantability or fitness for a particular purpose. No warranty may be created or extended by sales representatives or written sales material. The advice and strategies contained herein may not be suitable for your situation. Neither the Publisher nor author shall be liable for any loss of profit or any other commercial damages, including but not limited to special, incidental, consequential or other damages.

455 Boleskine Road, Victoria B.C. Canada, V8Z 1E7
info@ViVidroLabs.com | (250) 388-3531 | www.ViVidroLabs.com

Copyright © 2018 by ViVidro Labs Inc. All rights reserved.

 **ViVidro Labs Inc.**
First in Cardiovascular Device Testing