

PUMP HEAD

Product Codes: 16796

The ViVitro Pump Head can be used with the ViVitro SuperPump to create a circulatory, pulsatile flow loop. The Pump Head accommodates ventricular pressure monitoring using Luer fittings. It may be mounted directly to the SuperPump or can be used in conjunction with the ViVitro Viscoelastic Impedance Adaptor to attenuate pressure traces.

> Features

- Silicone ventricle membrane to isolate SuperPump piston from test fluid
- Low closing volume spring loaded disc valves
- 13 mm (1/2") hose connectors
- Easy to debubble
- Easy to disassemble and clean



CONTACT US

HEAD OFFICE

ViVitro Labs Inc.
455 Boleskine Road
Victoria, BC, V8Z 1E7, Canada
1 250.388.3531
info@vivitrolabs.com
skype: Vivitro.labs.inc.
www.vivitrolabs.com

NORTH AMERICA

Douglas Bigrigg, B. Eng.
North American Major Accounts Manager

ViVitro Labs Inc.
455 Boleskine Road
Victoria, BC, V8Z 1E7, Canada
1 250.388.3531 x217
DBigrigg@vivitrolabs.com
www.vivitrolabs.com

EUROPE, ISRAEL AND AFRICA

Jérémie Decressac
Business Development Manager

ViVitro Labs SASU
24-26 Boulevard Gay Lussac
Immeuble Energie
13014 Marseille, France
+33 4 86 68 68 10
JDecressac@vivitrolabs.com
www.vivitrolabs.com

CHINA

Shining Y. W. Sun
沈阳福帝斯科技有限公司

Shenyang Fortis Tech. Ltd.
024-23286966
yuwei.sun@fortis-hitech.com
www.fortis-hitech.com

INDIA

Ariel Lobo
Emergent Medi-Tech Ltd.
116-119 Damji Shamji Industrial Complex
9, L.B.S. Marg, Kurla (West)
Mumbai, 400070, India
91.22.67825555 x5
info@emergent.in
www.emergent.in

MIDDLE EAST

Mohan Das
Middle East Medical Supplies
PO Box 26143 Tahir Musbah Al Kendi Bldg.
Abu Dhabi, United Arab Emirates
+971 506221598
mohanda@memsuae.com
www.memsuae.com

No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form or by any means, electronic, mechanical, photocopying, recording, scanning or otherwise without the prior written permission of the Publisher.

Limit of Liability/Disclaimer of Warranty: While the publisher and author have used their best efforts in preparing this document, they make no representation or warranties with respect to the accuracy or completeness of the contents of this document and specifically disclaim any implied warranties of merchantability or fitness for a particular purpose. No warranty may be created or extended by sales representatives or written sales material. The advice and strategies contained herein may not be suitable for your situation. Neither the Publisher nor author shall be liable for any loss of profit or any other commercial damages, including but not limited to special, incidental, consequential or other damages.

455 Boleskine Road, Victoria B.C. Canada, V8Z 1E7 | (250) 388-3531
24-26 Blvd. Gay Lussac, Immeuble Energie, 13014 Marseille, France | +33 4 86 68 68 10
info@ViVitroLabs.com | www.ViVitroLabs.com

Copyright © 2021 by ViVitro Labs Inc. All rights reserved.

ViVitro
Labs