

# CUSTOMIZED CONDUIT MOUNT

Product Codes: 18903

*The ViVitro Conduit Mount (ACM) is used to mount various compliant conduits (native or synthetic root) to test surgical or percutaneous devices in the ViVitro Model Heart system. The ACM also enables valves to be sutured and mounted to native or synthetic roots in an anatomically correct fashion.*

Customized ACM accessories can be designed to fit device sizes and geometries. Custom designs can be created for items such as compliant conduits, as well as inflow and outflow hose barbs positioned on a semi-adjustable mount to accommodate various heights/length.



*Please contact ViVitro to discuss and customize an ACM for your unique application and testing requirements.*

No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form or by any means, electronic, mechanical, photocopying, recording, scanning or otherwise without the prior written permission of the Publisher.

**Limit of Liability/Disclaimer of Warranty:** While the publisher and author have used their best efforts in preparing this document, they make no representation or warranties with respect to the accuracy or completeness of the contents of this document and specifically disclaim any implied warranties of merchantability or fitness for a particular purpose. No warranty may be created or extended by sales representatives or written sales material. The advice and strategies contained herein may not be suitable for your situation. Neither the Publisher nor author shall be liable for any loss of profit or any other commercial damages, including but not limited to special, incidental, consequential or other damages.

455 Boleskine Road, Victoria B.C. Canada, V8Z 1E7 | (250) 388-3531  
24-26 Blvd. Gay Lussac, Immeuble Energie, 13014 Marseille, France | +33 4 86 68 68 10

info@ViVitroLabs.com | [www.ViVitroLabs.com](http://www.ViVitroLabs.com)

Copyright © 2022 by ViVitro Labs Inc. All rights reserved.

**ViVitro**  
Labs