

HEAT BATH

Product Code: 10679 N.A., 10680 Int'l

The Heat Bath Temperature Controller delivers rapid fluid heat-up, optimized for the ViVITRO Pulse Duplicator. These heated circulating baths feature a long life, high-wattage heater.

- Over Temperature and Low Reservoir level automatic system cut-offs eliminate risk of over heating system components

> Specifications

Reservoir Volume	6 litres (1.5 gallons)
Temperature Range	Ambient 15°C to 150°C
Stability	Digital $\pm 0.05^\circ\text{C}$
Work Area Dimensions	20 x 15 x 13 cm
Overall Dimensions	23 x 36 x 37 cm
Heater	1000 W
Pumping Capacity	12 L/min at 0 psig
Electrical Requirements	115 VAC, 60 Hz, 15 amp or 230 VAC, 50 Hz 10 amps (specify when ordering)



No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form or by any means, electronic, mechanical, photocopying, recording, scanning or otherwise without the prior written permission of the Publisher.

Limit of Liability/Disclaimer of Warranty: While the publisher and author have used their best efforts in preparing this document, they make no representation or warranties with respect to the accuracy or completeness of the contents of this document and specifically disclaim any implied warranties of merchantability or fitness for a particular purpose. No warranty may be created or extended by sales representatives or written sales material. The advice and strategies contained herein may not be suitable for your situation. Neither the Publisher nor author shall be liable for any loss of profit or any other commercial damages, including but not limited to special, incidental, consequential or other damages.

455 Boleskine Road, Victoria B.C. Canada, V8Z 1E7 | (250) 388-3531
24-26 Blvd. Gay Lussac, Immeuble Energie, 13014 Marseille, France | +33 4 86 68 68 10

info@ViVITRO Labs.com | www.ViVITRO Labs.com

Copyright © 2023 by ViVITRO Labs Inc. All rights reserved.

ViVITRO
Labs